


This guide will walk through the basic features of the My SmartLink Cloud to help get you started with using the program.

Dashboard

Upon first logging into the system, you will be presented with a Dashboard view of all patients in the system that may require intervention for various different needs including:

- Initial compliance risk at 3, 7, 30, and 45 days
- No data collected
- Therapy issues such as high AHI or leak
- Final compliance data for patient who met or did not meet current Medicare guidelines

Click “Review” to see specific patient details in calendar format for easy viewing.




Patient Information

From the main menu, you can view a patient by choosing “Patient’s”, then select the patient of choice and click “View Patient”

Inside the Patient Profile, you have access to view to following:


- View CPAP therapy reports including summary and nightly reports
- Write general or clinical notes in the patient file
- Add accessories such as masks, tubing, or water chambers
- Upload documents such as prescriptions or sleep lab reports

You can also add patient information such as general demographics information and contact details.



Physician Information:

Physicians can be added to the system and given access to view specific patients and make prescription changes if desired. Physicians are part of a global list to better enable physicians to see patients across multiple providers. To give physicians access to patient information, go into Patient Details in the Patient Profile and choose the Medical Care tab at the top



Data Management:

Data Management allows providers to upload SD cards into the My SmartLink Cloud. Once data is uploaded, the system will ask you to confirm the data to a specific patient and approve the upload.

